


## Quanser Interactive Labs Installation Guide

1. Go to <https://portal.quanser.com/Accounts/Login?returnUrl=/>
2. Register a new user using your McMaster email.



There is high chance that the confirmation email will be sent to spam. Please check your spam if can't see your confirmation email in your inbox.

3. Go to [portal.quanser.com/Downloads](https://portal.quanser.com/Downloads)
4. Download the **QLabs Virtual Experiments**. By clicking the “Download for Windows” Link. Note that this software only works with Windows OS.




5. After the download is complete, run the installer.
6. Test the software by launching the application



7. Select the Cube 2 – DC Motor



8. Select the servo Workspace



9. You should see a digital twin of the servomotor on your screen

